TECH BULLETIN

This is an addendum to all Aaron 1250 User's Guides and Service Guide. The revised memory feature and/or startup power and mode feature information is included in this bulletin:

EN

SECTION 1 USER'S GUIDE

- Memory

The unit automatically powers up to the Cut and Coag modes and their last selected power settings.

SECTION 4 USER'S GUIDE / SECTION 5 SERVICE GUIDE

Step 6

- The unit automatically powers up to the Cut and Coag modes and their last selected power settings. If the Blend, Fulguration or Bipolar mode is selected, the unit will default to that modes' last set power setting.

SECTION 4 USER'S GUIDE / SECTION 5 SERVICE GUIDE

The Aaron 1250U Electrosurgical Generator powers up to the Cut and Coag modes and their last set power settings. For example, you select the Cut mode and set the power to 50 watts. When you turn the unit off, it will default to Cut mode at 50 watts when its turned on again. The Coagulation mode will also default to the last set power setting. For the Blend, Fulguration or Bipolar modes and power settings, the unit will remember the last set power once the mode is again selected. For example, you select the Blend mode and set the power to 40 watts and select the Bipolar mode and set the power to 10 watts. When you turn the unit off, it will default to the Cut and Coag modes and their last selected power settings when its turned on again. If you then select Blend, the power will automatically default to 40 watts since it was your last set power setting.

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